

SEQUENCE LISTINGS

110> LIFECORD INC. et al.

<120> METHOD FOR ISOLATING AND CULTURING MULTIPOTENT PROGENITOR/STEM CELLS FROM UMBILICAL CORD BLOOD AND METHOD FOR INDUCING DIFFERENTIATION THEREOF

<130> PCA50101-LC1

<150> KR2004-6088

<151> 2004-01-30

<150> KR2005-6595

<151> 2005-01-25

<160> 28

<170> KopatentIn 1.71

<210> 1

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> forward primer specific for alkaline phosphatase

<400> 1

acgtggctaa gaatgcatc

20

<210> 2

<211> 19

<212> DNA

<213> Artificial Sequence

<220>

<223> reverse primer specific for alkaline phosphatase

<400> 2

2

ctggtaggcg atgtcctta

19

<210> 3

<211> 18

<212> DNA

<213> Artificial Sequence

<220>

<223> forward primer specific for type I procollagen

<400> 3

tgacgagacc aagaactg

18

<210> 4

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> reverse primer specific for type I procollagen

<400> 4

cgatccaaac cactgaaacc

20

<210> 5

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> forward primer specific for myoD

<400> 5

aatgtagcag gtgtaaccgt

20

<210> 6

3

<211> 20
<212> DNA
<213> Artificial Sequence

<220>
<223> reverse primer specific for myoD

<400> 6
gcctttattt tgatcacctg 20

<210> 7
<211> 20
<212> DNA
<213> Artificial Sequence

<220>
<223> forward primer specific for myogenin

<400> 7
cactacttct gtagcagggg 20

<210> 8
<211> 20
<212> DNA
<213> Artificial Sequence

<220>
<223> reverse primer specific for myogenin

<400> 8
tctctcaaac cgttcactt 20

<210> 9
<211> 20
<212> DNA
<213> Artificial Sequence

4

<220>

<223> forward primer specific for myosin heavy chain

<400> 9

tgtgaatgcc aaatgtgctt

20

<210> 10

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> reverse primer specific for myosin heavy chain

<400> 10

gtggagctgg gtatccttga

20

<210> 11

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> forward primer specific for Flt-1/VEGFR-1

<400> 11

ggtcttacgg agtattgctg

20

<210> 12

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> reverse primer specific for Flt-1/VEGFR-1

5

<400> 12
ctttcttttg ggtctctgtg 20

<210> 13
<211> 20
<212> DNA
<213> Artificial Sequence

<220>
<223> forward primer specific for KDR/VEGFR-2

<400> 13
ggacctggcg gcacgaaata 20

<210> 14
<211> 20
<212> DNA
<213> Artificial Sequence

<220>
<223> reverse primer specific for KDR/VEGFR-2

<400> 14
aggccggctc ttctgcttac 20

<210> 15
<211> 30
<212> DNA
<213> Artificial Sequence

<220>
<223> forward primer specific for ecNOS

<400> 15
aagacatttt cgggctcacg ctgcgcaccc 30

6

<210> 16
<211> 30
<212> DNA
<213> Artificial Sequence

<220>
<223> reverse primer specific for ecNOS

<400> 16
tggggtaggc actttagtag ttctcctaac 30

<210> 17
<211> 19
<212> DNA
<213> Artificial Sequence

<220>
<223> forward primer specific for VE-cadherin

<400> 17
gatgcagagg ctcatgatg 19

<210> 18
<211> 20
<212> DNA
<213> Artificial Sequence

<220>
<223> reverse primer specific for VE-cadherin

<400> 18
cttgcgactc acgcttgact 20

<210> 19
<211> 20
<212> DNA
<213> Artificial Sequence

<220>

<223> forward primer specific for vWF

<400> 19

caccgtttgc ccacccttcg

20

<210> 20

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> reverse primer specific for vWF

<400> 20

gccactggg agccgacact

20

<210> 21

<211> 19

<212> DNA

<213> Artificial Sequence

<220>

<223> forward primer specific for beta-actin

<400> 21

tgaaccaggc ttcagcatc

19

<210> 22

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> reverse primer specific for beta-actin

<400> 22
ggacttcgag caagatatgg 20

<210> 23
<211> 25
<212> DNA
<213> Artificial Sequence

<220>
<223> forward primer specific for HNF1-alpha

<400> 23
ttctaagctc agccagctgc agacg 25

<210> 24
<211> 25
<212> DNA
<213> Artificial Sequence

<220>
<223> reverse primer specific for HNF1-alpha

<400> 24
gctgaggttc tccggctctt tcaga 25

<210> 25
<211> 20
<212> DNA
<213> Artificial Sequence

<220>
<223> forward primer specific for cytokeratin-8

<400> 25
caatgccaag ctggaggatc 20

<210> 26
<211> 20
<212> DNA
<213> Artificial Sequence

<220>
<223> reverse primer specific for cytokeratin-8

<400> 26
acctcaggct ggcaatgact 20

<210> 27
<211> 24
<212> DNA
<213> Artificial Sequence

<220>
<223> forward primer specific for albumine

<400> 27
tgcttgaatg tgctgatgac aggg 24

<210> 28
<211> 25
<212> DNA
<213> Artificial Sequence

<220>
<223> reverse primer specific for albumine

<400> 28
aaggcaagtc agcaggcatc tcatc 25