

## SEQUENCE LISTINGS

110> LIFECORD INC. et al.

<120> METHOD FOR ISOLATING AND CULTURING MULTIPOTENT PROGENITOR/STEM CELLS FROM UMBILICAL CORD BLOOD AND METHOD FOR INDUCING DIFFERENTIATION THEREOF

<130> PCA50101-LC1

<150> KR2004-6088

<151> 2004-01-30

<150> KR2005-6595

<151> 2005-01-25

<160> 28

<170> KopatentIn 1.71

<210> 1

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> forward primer specific for alkaline phosphatase

<400> 1

acgtggcctaa gaatgtcatc

20

<210> 2

<211> 19

<212> DNA

<213> Artificial Sequence

<220>

<223> reverse primer specific for alkaline phosphatase

<400> 2

ctggtaggcg atgtcctta

19

<210> 3  
 <211> 18  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> forward primer specific for type I procollagen

<400> 3  
 tgacgagacc aagaactg

18

<210> 4  
 <211> 20  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> reverse primer specific for type I procollagen

<400> 4  
 cgatccaaac cactgaaacc

20

<210> 5  
 <211> 20  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> forward primer specific for myoD

<400> 5  
 aatgtagcag gtgtaaccgt

20

<210> 6

<211> 20  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> reverse primer specific for myoD

<400> 6  
gcctttatatt tgatcacctg 20

<210> 7  
<211> 20  
<212> DNA  
<213> Artificial Sequence  
  
<220>  
<223> forward primer specific for myogenin

<400> 7  
cactacttct gtagcagggg 20

<210> 8  
<211> 20  
<212> DNA  
<213> Artificial Sequence  
  
<220>  
<223> reverse primer specific for myogenin

<400> 8  
tctctcaaac cgtttcactt 20

<210> 9  
<211> 20  
<212> DNA  
<213> Artificial Sequence

<220>

<223> forward primer specific for myosin heavy chain

<400> 9

tgtgaatgcc aaatgtgctt

20

<210> 10

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> reverse primer specific for myosin heavy chain

<400> 10

gtggagctgg gtatccttga

20

<210> 11

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> forward primer specific for Flt-1/VEGFR-1

<400> 11

ggctttacgg agtattgctg

20

<210> 12

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> reverse primer specific for Flt-1/VEGFR-1

<400> 12  
ctttcttttg ggtctctgtg

20

<210> 13  
<211> 20  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> forward primer specific for KDR/VEGFR-2

<400> 13  
ggacctggcg gcacgaaata

20

<210> 14  
<211> 20  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> reverse primer specific for KDR/VEGFR-2

<400> 14  
aggccggctc tttcgcttac

20

<210> 15  
<211> 30  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> forward primer specific for ecNOS

<400> 15  
aagacatttt cgggctcacg ctgcgcaccc

30

<210> 16  
 <211> 30  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> reverse primer specific for ecNOS

<400> 16  
 tggggtaggc actttagtag ttctcctaac

30

<210> 17  
 <211> 19  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> forward primer specific for VE-cadherin

<400> 17  
 gatgcagagg ctcatgatg

19

<210> 18  
 <211> 20  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> reverse primer specific for VE-cadherin

<400> 18  
 ctgtcgactc acgcttgact

20

<210> 19  
 <211> 20  
 <212> DNA  
 <213> Artificial Sequence

<220>

<223> forward primer specific for vWF

<400> 19

caccgtttgc ccacccttcg

20

<210> 20

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> reverse primer specific for vWF

<400> 20

gccactggg agccgacact

20

<210> 21

<211> 19

<212> DNA

<213> Artificial Sequence

<220>

<223> forward primer specific for beta-actin

<400> 21

tgaaccaggc ttcagcatc

19

<210> 22

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> reverse primer specific for beta-actin

<400> 22  
ggacttcgag caagatatgg

20

<210> 23  
<211> 25  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> forward primer specific for HNF1-alpha

<400> 23  
ttctaagctc agccagctgc agacg

25

<210> 24  
<211> 25  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> reverse primer specific for HNF1-alpha

<400> 24  
gctgaggttc tccggctctt tcaga

25

<210> 25  
<211> 20  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> forward primer specific for cytokeratin-8

<400> 25  
caatgccaag ctggaggatc

20



<210> 26  
 <211> 20  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> reverse primer specific for cytokeratin-8

<400> 26  
 acctcaggct ggcaatgact

20

<210> 27  
 <211> 24  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> forward primer specific for albumine

<400> 27  
 tgcttgaatg tgctgatgac aggg

24

<210> 28  
 <211> 25  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> reverse primer specific for albumine

<400> 28  
 aaggcaagtc agcaggcatc tcatc

25

SEQUENCE LISTINGS

110> LIFECORD INC. et al.

<120> METHOD FOR ISOLATING AND CULTURING MULTIPOTENT PROGENITOR/STEM CELLS  
FROM UMBILICAL CORD BLOOD AND METHOD FOR INDUCING DIFFERENTIATION THEREOF

<130> Q96125

<150> KR2004-6088

<151> 2004-01-30

<150> KR2005-6595

<151> 2005-01-25

<150> PCT/KR2005/000278

<151> 2005-01-31

<160> 28

<170> KopatentIn 1.71

<210> 1

<211> 20

<212> DNA

<213> Artificial Sequence

<220>

<223> forward primer specific for alkaline phosphatase

<400> 1

acgtggctaa gaatgtcatc

20

<210> 2

<211> 19

<212> DNA

<213> Artificial Sequence

<220>

<223> reverse primer specific for alkaline phosphatase

<400> 2

ctggtaggcg atgtcctta

19

<210> 3

<211> 18

<212> DNA

<213> Artificial Sequence

<220>

<223> forward primer specific for type I procollagen

<400> 3

tgacgagacc aagaactg

18

<210> 4  
<211> 20  
<212> DNA  
<213> Artificial Sequence  
  
<220>  
<223> reverse primer specific for type I procollagen

<400> 4  
cgatccaaac cactgaaacc 20

<210> 5  
<211> 20  
<212> DNA  
<213> Artificial Sequence  
  
<220>  
<223> forward primer specific for myoD

<400> 5  
aatgtagcag gtgtaaccgt 20

<210> 6  
<211> 20  
<212> DNA  
<213> Artificial Sequence  
  
<220>  
<223> reverse primer specific for myoD

<400> 6  
gcctttattt tgatcacctg 20

<210> 7  
<211> 20  
<212> DNA  
<213> Artificial Sequence  
  
<220>  
<223> forward primer specific for myogenin

<400> 7  
cactacttct gtagcagggg 20

<210> 8  
<211> 20  
<212> DNA  
<213> Artificial Sequence  
  
<220>  
<223> reverse primer specific for myogenin

<400> 8  
tctctcaaac cgtttcactt 20

<210> 9  
<211> 20  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> forward primer specific for myosin heavy chain

<400> 9  
tgtgaatgcc aaatgtgctt 20

<210> 10  
<211> 20  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> reverse primer specific for myosin heavy chain

<400> 10  
gtggagctgg gtatccttga 20

<210> 11  
<211> 20  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> forward primer specific for Flt-1/VEGFR-1

<400> 11  
ggtcttacgg agtattgctg 20

<210> 12  
<211> 20  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> reverse primer specific for Flt-1/VEGFR-1

<400> 12  
ctttcttttg ggtctctgtg 20

<210> 13  
<211> 20

<212> \* DNA  
 <213> Artificial Sequence  
  
 <220>  
 <223> forward primer specific for KDR/VEGFR-2  
  
 <400> 13  
 ggacctggcg gcacgaaata 20  
  
 <210> 14  
 <211> 20  
 <212> DNA  
 <213> Artificial Sequence  
  
 <220>  
 <223> reverse primer specific for KDR/VEGFR-2  
  
 <400> 14  
 aggccggctc ttctgcttac 20  
  
 <210> 15  
 <211> 30  
 <212> DNA  
 <213> Artificial Sequence  
  
 <220>  
 <223> forward primer specific for ecNOS  
  
 <400> 15  
 aagacatttt cgggctcacg ctgcgcaccc 30  
  
 <210> 16  
 <211> 30  
 <212> DNA  
 <213> Artificial Sequence  
  
 <220>  
 <223> reverse primer specific for ecNOS  
  
 <400> 16  
 tggggtaggc acttttagtag ttctcctaac 30  
  
 <210> 17  
 <211> 19  
 <212> DNA  
 <213> Artificial Sequence  
  
 <220>  
 <223> forward primer specific for VE-cadherin  
  
 <400> 17

gatgca~~g~~agg ctcattgatg 19

<210> 18  
<211> 20  
<212> DNA  
<213> Artificial Sequence  
  
<220>  
<223> reverse primer specific for VE-cadherin

<400> 18  
cttgcgactc acgcttgact 20

<210> 19  
<211> 20  
<212> DNA  
<213> Artificial Sequence  
  
<220>  
<223> forward primer specific for vWF

<400> 19  
caccgtttgc ccacccttcg 20

<210> 20  
<211> 20  
<212> DNA  
<213> Artificial Sequence  
  
<220>  
<223> reverse primer specific for vWF

<400> 20  
gcccactggg agccgacact 20

<210> 21  
<211> 19  
<212> DNA  
<213> Artificial Sequence  
  
<220>  
<223> forward primer specific for beta-actin

<400> 21  
tgaaccaggc ttcagcatc 19

<210> 22  
<211> 20  
<212> DNA  
<213> Artificial Sequence

<220> '  
 <223> reverse primer specific for beta-actin

<400> 22  
 ggacttcgag caagatatgg 20

<210> 23  
 <211> 25  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> forward primer specific for HNF1-alpha

<400> 23  
 ttctaagctc agccagctgc agacg 25

<210> 24  
 <211> 25  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> reverse primer specific for HNF1-alpha

<400> 24  
 gctgagggttc tccggctctt tcaga 25

<210> 25  
 <211> 20  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> forward primer specific for cytokeratin-8

<400> 25  
 caatgccaag ctggaggatc 20

<210> 26  
 <211> 20  
 <212> DNA  
 <213> Artificial Sequence

<220>  
 <223> reverse primer specific for cytokeratin-8

<400> 26  
 acctcaggct ggcaatgact 20

<210> \* 27  
<211> 24  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> forward primer specific for albumine

<400> 27  
tgcttgaatg tgctgatgac aggg 24

<210> 28  
<211> 25  
<212> DNA  
<213> Artificial Sequence

<220>  
<223> reverse primer specific for albumine

<400> 28  
aaggcaagtc agcaggcatc tcac 25